



Liquid Biopsy Research - "We are where you start"

# HighPrep™ Pathogen Maxi Kit

Manual Revision v1.1

Catalog Nos. HPPM-D16, HPPM-D96, HPPM-D96X4

- Genomic DNA isolation from various pathogens from tissue, urine, serum, and fecal samples
- Magnetic beads based chemistry

## PROTOCOL

### Contents

Product Description and Process.....	1
Kit Contents, Storage, Stability.....	1
Preparation of Reagents .....	2
HighPrep Pathogen Maxi kit: Tissue Protocol.....	3
HighPrep Pathogen Maxi kit: Serum and Stool Protocol.....	6
HighPrep Pathogen Maxi kit: Urine Protocol.....	9
Troubleshooting Guide.....	12
Ordering and Related Product Information.....	13

#### **For Research Use Only. Not for use in diagnostic procedures.**

Information in this document is subject to change without notice.

MAGBIO GENOMICS, INC. DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. TO THE FULLEST EXTENT ALLOWED BY LAW, IN NO EVENT SHALL MAGBIO GENOMICS, INC. BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT MAGBIO GENOMICS, INC. IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

#### **TRADEMARKS**

The trademarks mentioned herein are the property of MagBio Genomics, Inc. or their respective owners.

## Product Description

MagBio Pathogen Maxi is paramagnetic bead-based kit, it is designed for efficient isolation of high-quality genomic DNA from various pathogens including gram positive and negative bacterial DNA, fungal spore DNA, and viral DNA and RNA from tissue, urine, serum, and fecal samples.

The kit's utilization method is adaptable on many major liquid handling workstations for high throughput processing, and will generate genomic DNA, suitable for PCR, restriction digestion, and hybridization applications.

## Process

The HighPrep™ Pathogen Maxi Kit uses a simple 3 step procedure: Lyse+Bind-Wash-Elute. Samples are lysed and DNA binds to the MAG-RQ1 magnetic beads in one step. Utilizing a magnetic separation device, the bound genomic DNA is separated from the solution and is washed. The final step is elution of high quality genomic DNA from the magnetic beads.

## Kit Contents and Storage

HighPrep™ Pathogen Maxi Kit Catalog No.	HPPM-D16	HPPM-D96	HPPM-D384	STORAGE
Number of Preps	16	96	384	
PL Buffer	10 ml	60 ml	240 ml	15-25°C
PB Buffer	1.5 ml	8 ml	30 ml	15-25°C
PP Buffer	4.5 ml	25 ml	100 ml	15-25°C
PBG Buffer	7 ml	40 ml	160 ml	15-25°C
PBB Buffer	7 ml	40 ml	160 ml	15-25°C
PBH Buffer	15 ml	88 ml	3 x 88 ml	15-25°C
MQ Wash Buffer	5 ml	30 ml	4 x 30 ml	15-25°C
MB Elution Buffer	2.5 ml	15 ml	50 ml	15-25°C
Pro K Solution <sup>2</sup>	380 µl	2.2 ml	9 ml	2-8°C
MAG-RQ1 Particles	380 µl	2.2 ml	9 ml	2-8°C

<sup>1</sup> Ethanol must be added prior to use. See Preparation of Reagents

## Stability

All components are stable for 12 months when stored accordingly.

<sup>2</sup>Pro K Solution comes in a ready to use solution. Component is stable for 1 year when stored at 15-25°C. For storage longer than 1 year, storage at 2-8°C is recommended.

## Safety Information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate material safety data sheets (MSDSs). MSDS can be downloaded from the "Product Resource" tab when viewing the product kit.

## Preparation of Reagents

Prepare the following components for each kit before use:

Catalog No.	Component	Add 100% Ethanol	Storage
<b>HPPM-D16</b>	MQ Wash Buffer	17.5 ml	Room Temp 15-25°C
<b>HPPM-D16</b>	PBH Buffer	28 ml	Room Temp 15-25°C

Components are stable for 1 year when stored closed at room temperature

Catalog No.	Component	Add 100% Ethanol	Storage
<b>HPPM-D96</b>	MQ Wash Buffer	70 ml	Room Temp 15-25°C
<b>HPPM-D96</b>	PBH Buffer	112 ml	Room Temp 15-25°C

Components are stable for 1 year when stored closed at room temperature

Catalog No.	Component	Add 100% Ethanol	Storage
<b>HPPM-D384</b>	MQ Wash Buffer	70 ml per Bottle	Room Temp 15-25°C
<b>HPPM-D384</b>	PBH Buffer	112 ml per Bottle	Room Temp 15-25°C

Components are stable for 1 year when stored closed at room temperature

## HighPrep™ Pathogen Maxi Kit : Tissue Protocol

***The following protocol can be applied for saliva samples.***

### Materials and Equipment to be Supplied by User

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate material safety data sheets (MSDSs) from each product supplier.


- Centrifuge capable of at least 3,500 x g
- Magnetic Separation Device (MBMS-10)
- Incubator capable of 65°C
- 1.5ml tubes
- Vortexer
- 100% ethanol
- Nuclease Free Water

### Things to do before starting

- Prepare PBH Buffer and MQ Wash Buffer according to the “Preparing Reagents” section on Page 2.
- Set an incubator to 65°C.
- Heat Elution Buffer to 70°C

### Protocol

1. **Add 25-30 mg tissue to each 1.5mL tube after mincing the tissue.**
2. **Add 525 µL PL Buffer to each sample.**
3. **Add 53µL PB Buffer and 20 µL Proteinase K Solution to each sample.**
4. **Vortex for 60 seconds to mix thoroughly.**
5. **Incubate at 65°C for 30 minutes. Mix once during incubation.**
6. **Centrifuge at 3,500 x g for 10 minutes.**
7. **Transfer 300 µL cleared supernatant to a 1.5mL Tube**
8. **Add 300 µL PBG Buffer, 300 µL PBB Buffer, and 20 µL Mag-RQ1 Particles to each sample. Vortex to mix thoroughly or pipet up and down 20 times.**

 Shake thoroughly the MAG-Bind particles to fully resuspend before use.

9. **Let sit at room temperature for 10 minutes.**

10. Place the 1.5mL tube on the Magnetic Separation Device to magnetize the Mag-RQ1 Particles. Let sit at room temperature until the Mag-RQ1 Particles are completely cleared from solution.
11. Aspirate and discard the cleared supernatant. Do not disturb the Mag-RQ1 Particles.
12. Remove the tube containing the Mag-RQ1 Particles from the Magnetic Separation Device.
13. Add 600 µL PBH Buffer to each sample. Resuspend the Mag-RQ1 Particles by vortexing or pipetting up and down 20 times.  
Note: PBH Buffer must be diluted with 100% ethanol prior to use. Please see Page 2 for instructions.
14. Let sit at room temperature for 2 minutes.
15. Place the 1.5mL tube on the Magnetic Separation Device to magnetize the Mag-RQ1 Particles. Let sit at room temperature until the Mag-RQ1 Particles are completely cleared from solution.
16. Aspirate and discard the cleared supernatant. Do not disturb the Mag-RQ1 Particles.
17. Remove the tube containing the Mag-RQ1 Particles from the Magnetic Separation Device.
18. Repeat Steps 13-17 once for a second VHB Wash step.
19. Add 600 µL MQ Wash Buffer to each sample. Resuspend the Mag-RQ1 Particles by vortexing or pipetting up and down 20 times.  
Note: MQ Wash Buffer must be diluted with 100% ethanol prior to use. Please see Page 2 for instructions.
20. Let sit at room temperature for 2 minutes.
21. Place the 1.5mL tube on the Magnetic Separation Device to magnetize the Mag-RQ1 Particles. Let sit at room temperature until the Mag-RQ1 Particles are completely cleared from solution.
22. Aspirate and discard the cleared supernatant. Do not disturb the Mag-RQ1 Particles.
23. Leave the tube on the Magnetic Separation Device.
24. Add 500 µL nuclease-free water (not provided) to each sample. Immediately aspirate the nuclease free water. Do not let the samples stay in contact with the nuclease-free water for more than 60 seconds.  
Note: If using an automated platform, use the maximum volume the tips will allow up to 600 µL.
25. Add 50-100 µL MB Elution Buffer heated to 70°C to each sample. Resuspend the Mag-RQ1 Particles by vortexing or pipetting up and down 20 times.
26. Let sit at room temperature for 5 minutes.
27. Place the 1.5mL tube on the Magnetic Separation Device to magnetize the Mag-RQ1 Particles. Let sit at room temperature until the Mag-RQ1 Particles are completely cleared

**from solution.**

- 28. Transfer the cleared supernatant containing purified DNA to a clean 1.5mL tube. Store the DNA at -20°C.**

## HighPrep™ Pathogen Maxi Kit : Serum and Stool Protocol

### Materials and Equipment to be Supplied by User


When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate material safety data sheets (MSDSs) from each product supplier.

- Centrifuge capable of at least 3,500 x g
- Magnetic Separation Device (MBMS-10)
- Incubator capable of 65°C
- 1.5ml tubes
- Vortexer
- 100% ethanol
- Nuclease Free Water

### Things to do before starting

- Prepare PBH Buffer and MQ Wash Buffer according to the “Preparing Reagents” section on Page 2.
- Set an incubator to 70°C.
- Heat Elution Buffer to 70°C

### Protocol

- 1. Add 250 µL serum or stool samples to each 1.5mL tube. If stool sample is solid, resuspend to 10% wt/volume in PBS before starting.**
- 2. Add 275 µL PL Buffer to each sample.**
- 3. Add 50µL PB Buffer and 20 µL Proteinase K Solution to each sample.**
- 4. Vortex for 60 seconds to mix thoroughly.**
- 5. Incubate at 65°C for 30 minutes. Mix once during incubation.**
- 6. Centrifuge at 3,500 x g for 10 minutes.**
- 7. Transfer 300 µL cleared supernatant to a 1.5mL Tube.**
- 8. Add 300 µL PBG Buffer, 300 µL PBB Buffer, and 20 µL Mag-RQ1 Particles to each sample. Vortex to mix thoroughly or pipet up and down 20 times.**  
  
 Shake thoroughly the MAG-Bind particles to fully resuspend before use.
- 9. Let sit at room temperature for 10 minutes.**
- 10. Place the 1.5mL tube on the Magnetic Separation Device to magnetize the Mag-RQ1**

**Particles. Let sit at room temperature until the Mag-RQ1 Particles are completely cleared from solution.**

- 11. Aspirate and discard the cleared supernatant. Do not disturb the Mag-RQ1 Particles.**
- 12. Remove the tube containing the Mag-RQ1 Particles from the Magnetic Separation Device.**
- 13. Add 600 µL PBH Buffer to each sample. Resuspend the Mag-RQ1 Particles by vortexing or pipetting up and down 20 times.**  
**Note: PBH Buffer must be diluted with 100% ethanol prior to use. Please see Page 2 for instructions.**
- 14. Let sit at room temperature for 2 minutes.**
- 15. Place the 1.5mL tube on the Magnetic Separation Device to magnetize the Mag-RQ1 Particles. Let sit at room temperature until the Mag-RQ1 Particles are completely cleared from solution.**
- 16. Aspirate and discard the cleared supernatant. Do not disturb the Mag-RQ1 Particles.**
- 17. Remove the tube containing the Mag-RQ1 Particles from the Magnetic Separation Device.**
- 18. Repeat Steps 13-17 once for a second PBH Wash step.**
- 19. Add 600 µL MQ Wash Buffer to each sample. Resuspend the Mag-RQ1 Particles by vortexing or pipetting up and down 20 times.**  
**Note: MQ Wash Buffer must be diluted with 100% ethanol prior to use. Please see Page 2 for instructions.**
- 20. Let sit at room temperature for 2 minutes.**
- 21. Place the 1.5mL tube on the Magnetic Separation Device to magnetize the Mag-RQ1 Particles. Let sit at room temperature until the Mag-RQ1 Particles are completely cleared from solution.**
- 22. Aspirate and discard the cleared supernatant. Do not disturb the Mag-RQ1 Particles.**
- 23. Leave the tube on the Magnetic Separation Device.**
- 24. Add 500 µL nuclease-free water (not provided) to each sample. Immediately aspirate the nuclease free water. Do not let the samples stay in contact with the nuclease-free water for more than 60 seconds.**  
**Note: If using an automated platform, use the maximum volume the tips will allow up to 600 µL.**
- 25. Add 50-100 µL MB Elution Buffer heated to 70°C to each sample. Resuspend the Mag-RQ1 Particles by vortexing or pipetting up and down 20 times.**
- 26. Let sit at room temperature for 5 minutes.**
- 27. Place the 1.5mL tube on the Magnetic Separation Device to magnetize the Mag-RQ1**



**Particles. Let sit at room temperature until the Mag-RQ1 Particles are completely cleared from solution.**

- 28. Transfer the cleared supernatant containing purified DNA to a clean 1.5mL tube. Store the DNA at -20°C.**

## HighPrep™ Pathogen Maxi Kit : Urine Protocol

### Materials and Equipment to be Supplied by User


When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate material safety data sheets (MSDSs) from each product supplier.

- Centrifuge capable of at least 3,500 x g
- Magnetic Separation Device (MBMS-10)
- Incubator capable of 65°C
- 1.5ml tubes
- Vortexer
- 100% ethanol
- Nuclease Free Water

### Things to do before starting

- Prepare PBH Buffer and MQ Wash Buffer according to the “Preparing Reagents” section on Page 2.
- Set an incubator to 70°C.
- Heat Elution Buffer to 70°C
- Prepare an Ice bucket

### Protocol

- 1. Add 250 µL urine sample to each 1.5mL tube.**
- 2. Add 275 µL PL Buffer to each sample.**
- 3. Add 50µL PB Buffer and 20 µL Proteinase K Solution to each sample.**
- 4. Vortex for 60 seconds to mix thoroughly.**
- 5. Incubate at 65°C for 30 minutes. Mix once during incubation.**
- 6. Add 200 µL PP Buffer to each well. Place the plate on ice for 5 minutes**
- 7. Centrifuge at 3,500 x g for 10 minutes.**
- 8. Transfer 300 µL cleared supernatant to a 1.5mL tube.**
- 9. Add 300 µL PBG Buffer, 300 µL RBB Buffer, and 20 µL Mag-RQ1 Particles to each sample. Vortex to mix thoroughly or pipet up and down 20 times.**
  -  Shake thoroughly the MAG-Bind particles to fully resuspend before use.
- 10. Let sit at room temperature for 10 minutes.**

11. Place the 1.5mL tube on the Magnetic Separation Device to magnetize the Mag-RQ1 Particles. Let sit at room temperature until the Mag-RQ1 Particles are completely cleared from solution.
12. Aspirate and discard the cleared supernatant. Do not disturb the Mag-RQ1 Particles.
13. Remove the tube containing the Mag-RQ1 Particles from the Magnetic Separation Device.
14. Add 600 µL PBH Buffer to each sample. Resuspend the Mag-RQ1 Particles by vortexing or pipetting up and down 20 times.  
Note: PBH Buffer must be diluted with 100% ethanol prior to use. Please see Page 2 for instructions.
15. Let sit at room temperature for 2 minutes.
16. Place the 1.5mL tube on the Magnetic Separation Device to magnetize the Mag-Bind® Particles RQ. Let sit at room temperature until the Mag-RQ1 Particles are completely cleared from solution.
17. Aspirate and discard the cleared supernatant. Do not disturb the Mag-RQ1 Particles.
18. Remove the tube containing the Mag-RQ1 Particles from the Magnetic Separation Device.
19. Repeat Steps 13-17 once for a second PBH Wash step.
20. Add 600 µL MQ Wash Buffer to each sample. Resuspend the Mag-RQ1 Particles by vortexing or pipetting up and down 20 times.  
Note: MQ Wash Buffer must be diluted with 100% ethanol prior to use. Please see Page 2 for instructions.
21. Let sit at room temperature for 2 minutes.
22. Place the 1.5mL tube on the Magnetic Separation Device to magnetize the Mag-RQ1 Particles. Let sit at room temperature until the Mag-RQ1 Particles are completely cleared from solution.
23. Aspirate and discard the cleared supernatant. Do not disturb the Mag-RQ1 Particles.
24. Leave the tube on the Magnetic Separation Device.
25. Add 500 µL nuclease-free water (not provided) to each sample. Immediately aspirate the nuclease free water. Do not let the samples stay in contact with the nuclease-free water for more than 60 seconds.  
Note: If using an automated platform, use the maximum volume the tips will allow up to 600 µL.
26. Add 50-100 µL MB Elution Buffer heated to 70°C to each sample. Resuspend the Mag-RQ1 Particles by vortexing or pipetting up and down 20 times.
27. Let sit at room temperature for 5 minutes.

- 28. Place the 1.5mL tube on the Magnetic Separation Device to magnetize the Mag-RQ1 Particles. Let sit at room temperature until the Mag-RQ1 Particles are completely cleared from solution.**
- 29. Transfer the cleared supernatant containing purified DNA to a clean 1.5mL tube. Store the DNA at -20°C.**

## Troubleshooting guide

Please use this guide to troubleshoot any problems that may arise. For further assistance, please contact technical support via:

Phone: 1-855-262-4246 (in US), outside US, 1-301-302-0144

Email: support@magbiogenomics.com

Symptoms	Possible Causes	Comments
A260/A230 ratio is low	Salt contamination	<ul style="list-style-type: none"> <li>Repeat the DNA isolation with a new sample.</li> <li>Extend the incubation time with VHB Buffer.</li> <li>Wash the Mag-Bind® Particles RQ with ethanol.</li> </ul>
A260/A280 ratio is high	RNA contamination	The protocol does not remove RNA. If desired, add 5 µL RNase A (25 mg/mL) after lysate is cleared and before binding buffers are added. Let sit at room temperature for 5 minutes.
Low DNA Yield or no DNA Yield	Poor homogenization of sample	Repeat the DNA isolation with a new sample, be sure to mix the sample with SLX-Mlus Buffer thoroughly. Use a commercial homogenizer if possible.
	DNA washed off	Make sure VHB Buffer and SPM Wash Buffer are mixed with ethanol.
	Water Wash extended	Make sure that water wash step does not exceed 60 seconds and the Mag-Bind® Particles RQ are not resuspended
	Mag-Bind® Particles RQ lost in process	After water is added during wash step Mag-Bind® Particles RQ will go into solution. Let magnetic beads remagnetize prior to aspirating liquid.
Problems in downstream applications	BSA not added to PCR mixture	Add BSA to a final concentration of 0.1 µg/mL to the PCR mixture
	Too much DNA inhibits PCR reactions	Dilute the DNA elute used in the downstream application if possible.
	Non-specific bands in downstream PCR	Use hot-start Taq polymerase mixture
	Inhibitory substance in the eluted DNA.	Check the A260/A230 ratio. Dilute the elute to 1:50 if necessary.

## Ordering and Related Product Information

### Post PCR and Next Gen library prep clean up system

Catalog No.	Product
AC-60005	HighPrep PCR (5 mL)
AC-60050	HighPrep PCR (50 mL)
AC-60500	HighPrep PCR (500 mL)

### BigDye Sanger Sequencing Cleanup

Catalog No.	Product
DT-70005	HighPrep DTR (5 mL)
DT-70050	HighPrep DTR (50 mL)
DT-70500	HighPrep DTR (500 mL)

### Magnetic Separation Devices

Catalog No.	Description
MYMAG-96	Handheld Magnetic Separation Device (96 well microplate format)
MBMS-10	MagStip magnetic stand (1.5 mL x 10)
MBMS-31550	15ml and 50ml magnetic stand combo. (3x15ml and 3x50ml)

### cfDNA Purification Kit

Catalog No.	Product	Description	Preps
CFK-D10-400UL	CF-Kapture 21 Kit (200-400µl) 10 preps	Purification of cell-free DNA (cfDNA) from 200-400 µl STABILIZED plasma	10
CFK-D5-5ML	CF-Kapture 21 Kit (3-5ml) 5 preps	Purification of cell-free DNA (cfDNA) from 3-5 ml STABILIZED plasma	5
CFK-D50-400UL	CF-Kapture 21 Kit (200-400µl) (50 preps)	Purification of cell-free DNA (cfDNA) from 200-400 µl STABILIZED plasma	50
CFK-D50-2ML	CF-Kapture 21 Kit (1-2ml) 50 preps	Purification of cell-free DNA (cfDNA) from 1-2 ml STABILIZED plasma	50
CFK-D50-5ML	CF-Kapture 21 Kit (3-5 ml) 50 preps	Purification of cell-free DNA (cfDNA) from 3-5 ml STABILIZED plasma	50

### Whole blood stabilization tubes

Catalog No.	Product	Description
BS21-CF10-100	Blood STASIS 21-ccfDNA 9 mL (100)	100 tubes: 2 ml Additive, 7 ml blood draw volume
BS21-CF6-100	Blood STASIS 21-ccfDNA 6 mL (100)	100 tubes: 1.5 ml Additive, 4.5 ml blood draw volume
BS21-CF3-200	Blood STASIS 21-ccfDNA 3 mL (200)	200 tubes: 0.5 ml Additive, 2.5 ml blood draw volume

**RNA or cDNA for in vitro applications clean up system**

<b>Catalog No.</b>	<b>Product</b>
RC-60005	HighPrep RNA Elite (5 mL)
RC-60050	HighPrep RNA Elite (50 mL)
RC-60500	HighPrep RNA Elite (500 mL)













Liquid Biopsy Research - "We are where you start"

[www.magbiogenomics.com](http://www.magbiogenomics.com)