



# HighPrep PCR vs Leading SPRI Bead-Based Reagent

## DNA Clean-up and Recovery Rate




**Figure 1: Removal of 50 bp DNA fragment.** About 10  $\mu$ L of 50 bp DNA ladder was purified with HighPrep PCR and the Leading SPRI Bead-Based Reagent using a bead ratio of 1.8x following the manufacturer's protocol. Fragments above 100 bp were retained by both reagents and the 50 bp DNA fragment was removed.



**Figure 2: Side-by-side comparison of DNA recovery post cleanup using HighPrep PCR and the Leading SPRI Bead-Based Reagent.** About 10  $\mu$ L of 50 bp DNA ladder was purified using a bead ratio of 1.8x and DNA recovery was assessed. Both reagents had similar recovery.

## DNA Size Selection

Selection of 250 bp – 350 bp DNA fragment



**Figure 3: Efficient double-sided DNA size selection using HighPrep PCR as a direct replacement for the Leading SPRI Bead-Based Reagent.** Agilent TapeStation fragment analysis of 250-350 bp DNA size selection (0.75x/0.2x) from sheared genomic DNA. The same size selection ratio was applied for both reagents.

## Conclusion

- HighPrep PCR and the Leading SPRI Bead-Based Reagent demonstrate equivalent purification capabilities, ensuring effective DNA fragment retention while efficiently removing <50 bp fragments.
- HighPrep PCR delivers identical size selection as the Leading SPRI Bead-Based Reagent, making it a direct and cost-effective replacement.
- Stability assessments confirm that HighPrep PCR maintains efficiency over time, proving its reliability for long-term applications